

STANDARD BLOOD COLLECTION SET

STERILE, SINGLE-USE FOR BLOOD COLLECTION

PRODUCT DATASHEET



STANDARD BLOOD COLLECTION SET

STERILE, SINGLE-USE FOR BLOOD COLLECTION

- Ultra sharp three-beveled needle to minimize discomfort
- Siliconization of needle to minimize penetration and gliding force
- 20 G luer adapter with self-sheathing needle system
- Large, soft wings for easy grip
- Not made with latex, DEHP or PVC
- Available with 19 mm tubing length
- Available in 21 G, 22 G, 23 G
- For use by healthcare professionals only

PRODUCT COMPLIANCE

- Complies with the following norms, directives, and regulations:
 - EN ISO 14971
 - ISO 594-2
 - EN ISO 11607-2
 - ISO 11135-1
 - ISO 594-1
 - EN ISO 11607-1
 - ISO 9626
- Labels contain 12 languages

PRODUCT RANGE OVERVIEW

Reference Code	Description	Tubing length	Flow rate not less than	Priming value	Color
BCS-21-19-ET-FL	Blood Collection Set 21G x 3/4" (0.80 x 19 mm)	190 mm	≥ 15.0 ml/min	0.39 ml	Green
BCS-22-19-ET-FL	Blood Collection Set 22G x 3/4" (0.70 x 19 mm)	190 mm	≥ 12.0 ml/min	0.39 ml	Black
BCS-23-19-ET-FL	Blood Collection Set 23G x 3/4" (0.60 x 19 mm)	190 mm	≥ 7.0 ml/min	0.39 ml	Blue



MATERIALS USED

Lubricant

Material: Silicone

Tip design

Three-beveled

Needle (cannula)

Material: Stainless steel SUS-304

Wing hub

Material: Polyvinyl chloride (PVC)

Glue

Material: Epoxy (resin)

Rubber sleeve


Material: Synthetic rubber

Luer connector

Material: Acrylonitrile Butadiene
Styrene (ABS)

Tubing

Material: Polyvinyl chloride (PVC)



NIPRO ADDRESS

NIPRO CORPORATION (Headquarters)

3-9-3, Honjo-Nishi, Kita-ku, Osaka,
531-8510, Japan
Tel: +81 6 6372 2331 Fax: +81 6 6371 7422
<https://www.nipro.com.jp/en/>

NIPRO ASIA PTE LTD

6 Temasek Boulevard
#31-01/02 Suntec Tower Four
Singapore 038986
Tel: +65 6588 2281 Fax: +65 6588 2481
Email: niproasia@nipro.com.sg
<http://www.niproasia.com.sg/>

NIPRO MALAYSIA SDN. BHD

B-3-2 The Ascent Paradigm,
No.1 Jalan SS7/26A, Kelana Jaya,
47301 Petaling Jaya, Selangor, Malaysia
Tel: +603 7622 1190/1191/1192
Fax: +603 7622 180
Email: office@nipro.com.my

NIPRO JMI MEDICAL LTD

Unique Heights, Level 11,
117, Kazi Nazrul Islam Avenue, Ramna,
Dhaka-1217, Bangladesh
Tel: +880 2 55138723 Fax: +880 2 55138725
Email: info@jmigroup-bd.com

NIPRO AUSTRALIA PTY LTD

Suite 2. 02, 657 Pacific Highway,
St. Leonards, NSW 2065, Australia
Tel: +61 2 9437 9000 Fax: + 61 2 9437 0860
Email: support@niproaustralia.com.au
<http://www.niproaustralia.com.au/>

NIPRO MEDICAL PVT LTD

24-CCA, DHA Phase 8, (Ex Park View),
Lahore-54792, Pakistan
Tel: +92 4237175802/4237175803 Fax: +92 4237175804
Email: info@nipromedical.com.pk

NIPRO SALES (THAILAND) CO., LTD.

89/170 Nipro Building, Viphavadee-Rangsit Road,
Talad-Bangkhen, Laksi, Bangkok, THAILAND 10210
Tel : +66 2551 2000 Fax : +66 2551 2006
Email : sales@niprosales.com

NIPRO SALES VIETNAM COMPANY LIMITED

Floor 15, Unit 15.06 & 15.07, Viettel Complex Build-ing 285
Cach Mang Thang 8, Ward 12, District 10, Ho Chi Minh city,
Vietnam
Tel : +84 908 260 405
Email : nhuy@niprosales.com

NIPRO SALES (THAILAND) CO., LTD. -MYANMAR REPRESENTATIVE OFFICE

No. E 302, Manawhary Housing, Baho Road, Ahlone
Township, Yangon, Myanmar
Tel : +66 2551 2000 Fax : +66 2551 2006
Email : win@niprosales.com

NIPRO MEDICAL PHILIPPINES CORPORATION

19/F Unit C, Menarco Tower, 32nd Street, BGC, Taguig,
Metro Manila, Philippines
Tel : +63 285419362
Email : info@nipro.com.ph