

# HYPODERMIC NEEDLES

FOR INJECTION THERAPY

PRODUCT DATA SHEET



# HYPODERMIC NEEDLES

FOR INJECTION THERAPY

- Ultra sharp three-beveled needle
- Siliconization on lancet point minimizes penetration and gliding force
- Available in thin or regular wall
- Color-coded to conform to ISO 6009
- Stainless steel cannula protected by needle cap
- Not made with latex or PVC
- For use by healthcare professionals only

## PRODUCT COMPLIANCE

- |               |                    |               |
|---------------|--------------------|---------------|
| — ISO 594-1   | — ISO 2859-1       | — ISO 13485   |
| — ISO 594-2   | — ISO 6009         | — ISO 15223-1 |
| — EN 868-5    | — ISO 7864         |               |
| — CSN EN 1707 | — EN ISO 9626 + A1 |               |

## MATERIALS USED

### Cutting angle

Three-beveled

### Needle (cannula)

Material: Stainless steel SUS-304

### Siliconization of cannula

Material: Silicone

### Glue

Material: Epoxy (resin)

### Needle hub and needle cap

Material: Polypropylene (PP)



## PRODUCT RANGE OVERVIEW

Reference code	Description	Gauge (G)	Length (mm)	Color
HN-1625-ET	Hypo Ndl 16 G x 1" (1.6x25 mm) LR TW	16	25	
HN-1638-ET	Hypo Ndl 16 G x 1 1/2" (1.6x38 mm) LR TW	16	38	
HN-1825-ET	Hypo Ndl 18 G x 1" (1.2x25 mm) LR TW	18	25	
HN-1832-ET	Hypo Ndl 18 G x 1 1/4" (1.2x32 mm) LR TW	18	32	
HN-1838-ET	Hypo Ndl 18 G x 1 1/2" (1.2x38 mm) LR TW	18	38	
HN-1850-ET	Hypo Ndl 18 G x 2 (1.2x50 mm) LR TW	18	50	
HN-1925-ET	Hypo Ndl 19 G x 1" (1.1x25 mm) LR TW	19	25	
HN-1932-ET	Hypo Ndl 19 G x 1 1/4" (1.1x32 mm) LR TW	19	32	
HN-1938-ET	Hypo Ndl 19 G x 1 1/2" (1.1x38 mm) LR TW	19	38	
HN-1950-ET	Hypo Ndl. 19 G x 2" (1.1x50 mm) LR TW	19	50	
HN-2025-ET	Hypo Ndl 20 G x 1" (0.9x25 mm) LR TW	20	25	
HN-2038-ET	Hypo Ndl 20 G x 1 1/2" (0.9x38 mm) LR TW	20	28	
HN-2050-ET	Hypo Ndl. 20 G x 2" (0.9x50 mm) LR TW	20	50	
HN-2070-ET	Hypo Ndl. 20 G x 2 3/4" (0.9x70 mm) LR TW	20	70	
HN-2116-ET	Hypo Ndl 21 G x 5/8" (0.8x16 mm) LR TW	21	16	
HN-2125-ET	Hypo Ndl 21 G x 1" (0.8x25 mm) LR TW	21	25	
HN-2132-ET	Hypo Ndl 21 G x 1 1/4" (0.8x32 mm) LR TW	21	32	
HN-2138-ET	Hypo Ndl 21 G x 1 1/2" (0.8x38 mm) LR TW	21	38	
HN-2150-ET	Hypo Ndl. 21 G x 2" (0.8x50 mm) LR TW	21	50	
HN-2170-ET	Hypo Ndl. 21 G x 2 3/4" (0.8x70 mm) LR TW	21	70	
HN-2216-ET	Hypo Ndl 22 G x 5/8" (0.7x16 mm) LR TW	22	16	
HN-2219-ET	Hypo Ndl 22 G x 3/4" (0.7x19 mm) LR TW	22	19	
HN-2225-ET	Hypo Ndl 22 G x 1" (0.7x25 mm) LR TW	22	25	
HN-2232-ET	Hypo Ndl 22 G x 1 1/4" (0.7x32 mm) LR TW	22	32	
HN-2238-ET	Hypo Ndl 22 G x 1 1/2" (0.7x38 mm) LR TW	22	38	
HN-2250-ET	Hypo Ndl. 22 G x 2" (0.7x50 mm) LR TW	22	50	
HN-2316-ET	Hypo Ndl 23 G x 5/8" (0.6x16 mm) LR TW	23	16	
HN-2325-ET	Hypo Ndl 23 G x 1" (0.6x25 mm) LR TW	23	25	
HN-2332-ET	Hypo Ndl 23 G x 1 1/4" (0.6x32 mm) LR TW	23	32	
HN-2338-ET	Hypo Ndl 23 G x 1 1/2" (0.6x38 mm) LR TW	23	38	
HN-2425-ET	Hypo Ndl 24 G x 1" (0.55x25 mm) LR TW	24	25	
HN-2432-ET	Hypo Ndl 24 G x 1 1/4" (0.55x32 mm) LR TW	24	32	
HN-2438-ET	Hypo Ndl 24 G x 1 1/2" (0.55x38 mm) LR TW	24	38	
HN-2513-ET	Hypo Ndl 25 G x 1/2" (0.5x13 mm) LR TW	25	13	
HN-2516-ET	Hypo Ndl 25 G x 5/8" (0.5x16 mm) LR TW	25	16	
HN-2525-ET	Hypo Ndl 25 G x 1" (0.5x25 mm) LR TW	25	25	
HN-2532-ET	Hypo Ndl 25 G x 1 1/4" (0.5x32 mm) LR TW	25	32	
HN-2538-ET	Hypo Ndl 25 G x 1 1/2" (0.5x38 mm) LR TW	25	38	
HN-2613-ET	Hypo Ndl 26 G x 1/2" (0.45x13 mm) LR RW	26	13	
HN-2622-ET	Hypo Ndl 26 G x 7/8" (0.45x22 mm) LR RW	26	22	
HN-2713-ET	Hypo Ndl 27 G x 1/2" (0.4x13 mm) LR RW	27	13	
HN-2719-ET	Hypo Ndl 27 G x 3/4" (0.4x19 mm) LR RW	27	19	
HN-2732-ET	Hypo Ndl 27 G x 1 1/4" (0.4x32 mm) LR RW	27	32	
HN-2738-ET	Hypo Ndl 27 G x 1 1/2" (0.4x38 mm) LR RW	27	38	
HN-2913-ET	Hypo Ndl 29 G x 1/2" (0.34x13 mm) LR RW	29	13	
HN-3013-ET	Hypo Ndl 30 G x 1/2" (0.3x13 mm) LR TW	30	13	

## NIPRO ADDRESS

### **NIPRO CORPORATION (Headquarters)**

3-9-3, Honjo-Nishi, Kita-ku, Osaka,  
531-8510, Japan  
Tel: +81 6 6372 2331 Fax: +81 6 6371 7422  
<https://www.nipro.com.jp/en/>

### **NIPRO ASIA PTE LTD**

6 Temasek Boulevard  
#31-01/02 Suntec Tower Four  
Singapore 038986  
Tel: +65 6588 2281 Fax: +65 6588 2481  
Email: [niproasia@nipro.com.sg](mailto:niproasia@nipro.com.sg)  
<http://www.niproasia.com.sg/>

### **NIPRO MALAYSIA SDN. BHD**

B-3-2 The Ascent Paradigm,  
No.1 Jalan SS7/26A, Kelana Jaya,  
47301 Petaling Jaya, Selangor, Malaysia  
Tel: +603 7622 1190/1191/1192  
Fax: +603 7622 180  
Email: [office@nipro.com.my](mailto:office@nipro.com.my)

### **NIPRO JMI MEDICAL LTD**

Unique Heights, Level 11,  
117, Kazi Nazrul Islam Avenue, Ramna,  
Dhaka-1217, Bangladesh  
Tel: +880 2 55138723 Fax: +880 2 55138725  
Email: [info@jmigroup-bd.com](mailto:info@jmigroup-bd.com)

### **NIPRO AUSTRALIA PTY LTD**

Suite 2. 02, 657 Pacific Highway,  
St. Leonards, NSW 2065, Australia  
Tel: +61 2 9437 9000 Fax: + 61 2 9437 0860  
Email: [support@niproaustralia.com.au](mailto:support@niproaustralia.com.au)  
<http://www.niproaustralia.com.au/>

### **NIPRO MEDICAL PVT LTD**

24-CCA, DHA Phase 8, (Ex Park View),  
Lahore-54792, Pakistan  
Tel: +92 4237175802/4237175803 Fax: +92 4237175804  
Email: [info@nipromedical.com.pk](mailto:info@nipromedical.com.pk)

### **NIPRO SALES (THAILAND) CO., LTD.**

89/170 Nipro Building, Viphavadee-Rangsit Road,  
Talad-Bangkhen, Laksi, Bangkok, THAILAND 10210  
Tel : +66 2551 2000 Fax : +66 2551 2006  
Email : [sales@niprosales.com](mailto:sales@niprosales.com)

### **NIPRO SALES VIETNAM COMPANY LIMITED**

Floor 15, Unit 15.06 & 15.07, Viettel Complex Build-ing 285  
Cach Mang Thang 8, Ward 12, District 10, Ho Chi Minh city,  
Vietnam  
Tel : +84 908 260 405  
Email : [nhuy@niprosales.com](mailto:nhuy@niprosales.com)

### **NIPRO SALES (THAILAND) CO., LTD. -MYANMAR REPRESENTATIVE OFFICE**

No. E 302, Manawhary Housing, Baho Road, Ahlone  
Township, Yangon, Myanmar  
Tel : +66 2551 2000 Fax : +66 2551 2006  
Email : [win@niprosales.com](mailto:win@niprosales.com)

### **NIPRO MEDICAL PHILIPPINES CORPORATION**

19/F Unit C, Menarco Tower, 32nd Street, BGC, Taguig,  
Metro Manila, Philippines  
Tel : +63 285419362  
Email : [info@nipro.com.ph](mailto:info@nipro.com.ph)