

SYRINGES

- Smooth, accurate motion along the entire length of the barrel
- Clear barrel & precise volume reading
- Large finger grips for secure handling



Product Range

| Syringe Capacity | Lure Tip | | | Individual Packing | Package Inner/Carton |
|------------------|----------|------|----------|--------------------|----------------------|
| | Slip | Lock | Catheter | | |
| 1mL | ✓ | ✓ | | Wrapper or Blister | 100 pcs / 1000 pcs |
| 1mL Tuberculin | ✓ | | | Wrapper or Blister | 100 pcs / 1000 pcs |
| 3mL | ✓ | ✓ | | Wrapper or Blister | 100 pcs / 1000 pcs |
| 5mL | ✓ | ✓ | | Wrapper or Blister | 100 pcs / 1000 pcs |
| 10mL | ✓ | ✓ | | Wrapper or Blister | 100 pcs / 1000 pcs |
| 20mL | ✓ | ✓ | | Wrapper or Blister | 50 pcs / 500 pcs |
| 30mL | ✓ | ✓ | ✓ | Blister | 50 pcs / 400 pcs |
| 50mL | ✓ | ✓ | ✓ | Blister | 25 pcs / 200 pcs |

Product Range with Needle

| O.D Gauge | 20G | 21G | 22G | 23G | 24G | 25G | 26G | 27G |
|--------------------|-------------------------------------|------------|-------|-----------|---------------|--------|--------|-------------|
| O.D Metric | 0.9mm | 0.8mm | 0.7mm | 0.6mm | 0.55mm | 0.5mm | 0.45mm | 0.4mm |
| ISO Standard Color | Yellow | Deep Green | Black | Deep Blue | Medium Purple | Orange | Brown | Medium Grey |
| Length | 1/2", 5/8", 3/4", 1", 1-1/4, 1-1/2" | | | | | | | |



NIPRO CORPORATION (Headquarters)

3-9-3, Honjo-Nishi, Kita-ku, Osaka,
531-8510, Japan
Tel: +81 6 6372 2331 Fax: +81 6 6371 7422
<https://www.nipro.com.jp/en/>

NIPRO ASIA PTE LTD

6 Temasek Boulevard
#31-01/02 Suntec Tower Four
Singapore 038986
Tel: +65 6588 2281 Fax: +65 6588 2481
Email: niproasia@nipro.com.sg
<http://www.niproasia.com.sg/>

NIPRO MALAYSIA SDN. BHD

B-3-2 The Ascent Paradigm,
No.1 Jalan SS7/26A, Kelana Jaya,
47301 Petaling Jaya, Selangor, Malaysia
Tel: +603 7622 1190/1191/1192
Fax: +603 7622 180
Email: office@nipro.com.my

NIPRO JMI MEDICAL LTD

Unique Heights, Level 11,
117, Kazi Nazrul Islam Avenue, Ramna,
Dhaka-1217, Bangladesh
Tel: +880 2 55138723 Fax: +880 2 55138725
Email: info@jmgroupp-bd.com

NIPRO AUSTRALIA PTY LTD

Suite 2. 02, 657 Pacific Highway,
St. Leonards, NSW 2065, Australia
Tel: +61 2 9437 9000 Fax: + 61 2 9437 0860
Email: support@niproaustralia.com.au
<http://www.niproaustralia.com.au/>

NIPRO MEDICAL PVT LTD

24-CCA, DHA Phase 8, (Ex Park View),
Lahore-54792, Pakistan
Tel: +92 4237175802/4237175803 Fax: +92
4237175804
Email: info@nipromedical.com.pk

NIPRO SALES (THAILAND) CO., LTD.

89/170 Nipro Building, Viphavadee-Rangsit Road,
Talad-Bangkhen, Laksi, Bangkok, THAILAND 10210
Tel : +66 2551 2000 Fax : +66 2551 2006
Email : sales@niprosales.com

NIPRO SALES VIETNAM COMPANY LIMITED

Floor 15, Unit 15.06 & 15.07, Viettel Complex Build-
ing 285 Cach Mang Thang 8, Ward 12, District 10,
Ho Chi Minh city, Vietnam
Tel : +84 908 260 405
Email : nhuy@niprosales.com

**NIPRO SALES (THAILAND) CO., LTD. -
MYANMAR REPRESENTATIVE OFFICE**

No. E 302, Manawhary Housing, Baho Road, Ahlone
Township, Yangon, Myanmar
Tel : +66 2551 2000 Fax : +66 2551 2006
Email : win@niprosales.com