



Guide
Extension
Catheter

GUIDE PLUS™ II



Features



Excellent crossability

with braided tube & hydrophilic coating

The braided tubing enhances pushability while the hydrophilic coating realizes excellent deliverability.




Smooth device insertion

A unique shaped port provides for easy device insertion.

Catheter design


GUIDEPLUS II ST

Inner diameter



GUIDEPLUS II EL

Inner diameter



Ordering information

Product Name	Model No.	Maximum External Diameter	Distal Shaft Length	Compatible Guiding Catheter	Inner Diameter
GUIDEPLUS II	GUIDEPLUS 2 ST INT	1.76mm / 5.3F	250mm	6F / ≥ 0.070 inch	1.33mm / 0.052inch
	GUIDEPLUS 2 EL INT				1.43mm / 0.056inch



VASCULAR DIV.

http://med.nipro.co.jp/en_products



NIPRO CORPORATION
3-9-3, Honjo-Nishi, Kita-Ku, Osaka 531-8510, Japan

- ▶ Please consult your sales representative regarding stock availability.
- ▶ Please refer to the Instructions For Use for further details regarding product use and specifications.
- ▶ GUIDEPLUS is a trademark of NIPRO Corporation.

[Manufacturing facility]
NIPRO CORPORATION ODATE FACTORY
8-7, Hanukiyachi, Niida, Odate-shi, Akita, 018-5794 Japan

GPLUS II-2006-E7