

Micro
Catheter




Supporting
Your success

Features

Beyond simply reducing procedural burden by facilitating easy guidewire exchange, the low profile of Mogul ensures deliverability when needed to provide guidewire back up support for those difficult cases, and the catheter body internal design yields pushability required for lesion penetration of CTO's.

Catheter Design



Ordering information

Model No.	Outer Diameter		Compatible Guidewire	Radiopaque Marker	Working Shaft Length
	Distal	Proximal			
AMOG130S	0.58 mm (1.7F)	0.82 mm (2.5F)	0.014 inch (0.36 mm)	1	130 cm