

---


Guiding  
Catheter

---

# *Profit*<sup>™</sup> 6F Plus

Guiding Catheter

*Gain the advantage*



# Profit™ 6F Plus

Guiding Catheter

## Features

### 0.072 inch Inner Diameter

A large internal diameter (0.072 inch) and PTFE liner provides a low frictional coefficient, making it ideal for large devices and performing KBT technique.

### Superior control

While offering strong kink resistance and remarkable operability, the flexible shaft provides for required support during the procedure.






### Optimal backup force




A gradual transition in rigidity from tip to shaft ensures gentle ostial engagement while maintaining strong back-up force.

### Safe soft tip

Along with a flexible and radiopaque tip, a 'rounding' process increases the safety profile when engaging the vessel.

## Ordering information

Product Type	Shape	Model Number
<b>Judkins Right</b> 	JR	P6P-JR35
		P6P-JR40
		P6P-JR35SH
		P6P-JR40SH
<b>Radial Unlimited</b> 	RU	P6P-RU
		P6P-RUM375
		P6P-RUSH
<b>Radial Right</b> 	RR	P6P-RR
<b>Short Amplatz Left</b> 	SAL	P6P-SAL1
		P6P-SAL1SH
<b>Amplatz Left</b> 	AL	P6P-AL075
		P6P-AL1
		P6P-AL15
		P6P-AL2
		P6P-AL075SH
		P6P-AL1SH

Product Type	Shape	Model Number
<b>Judkins Left Short Tip</b> 	JL	P6P-JL35ST
		P6P-JL40ST
<b>Judkins Left</b> 	JL	P6P-JL35
		P6P-JL40
		P6P-JL45
<b>Special Support</b> 	SS	P6P-JL35SH
		P6P-JL40SH
		P6P-SS325
		P6P-SS35
		P6P-SS375
		P6P-SS40
		P6P-SS325SH
		P6P-SS35SH