

ACID CONCENTRATES

Bags, Canisters, Octabins



ACID CONCENTRATES

Bags, Canisters, Octabins

Ready-to-use liquid acid concentrates intended for use with powder bicarbonate concentrate.

Acid concentrate is an essential component in the preparation of dialysis fluid, whereby the hemodialysis machine mixes liquid acid concentrate with powder bicarbonate concentrate and purified water.

The standard mixing ratio of Nipro acid concentrates is 1:44.



PORTFOLIO

	Na ⁺	Ca ⁺⁺	Mg ⁺⁺	K ⁺	Cl ⁻	CH ₃ COO ⁻	HCO ₃ ⁻	Glucose	Osmolarity
	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	g/l	mOsm/l
306-A [A+B]	139	1.25	0.48	1.50	103.14	2.85	37.96	1.50	294.50
307-A [A+B]	139	1.75	0.48	1.50	104.14	2.85	37.96	1.50	296.00
308-A [A+B]	139	1.50	0.48	1.50	103.64	2.85	37.96	1.00	292.48
310-A [A+B]	139	1.50	0.50	1.00	103.20	2.85	37.96	1.00	291.55
313-A [A+B]	140	1.50	0.50	1.50	108.50	3.00	34.00	1.00	294.55
314-A [A+B]	140	1.25	0.50	1.50	108.00	3.00	34.00	1.00	293.80
316-A [A+B]	140	0.63	0.25	1.50	102.44	4.00	36.81	1.00	291.18
318-A [A+B]	140	1.75	0.50	2.00	109.50	3.00	34.00	1.50	299.08
319-A [A+B]	139	1.75	0.48	3.00	105.64	2.85	37.96	1.50	299.00
320-A [A+B]	139	1.25	0.50	3.00	104.70	2.85	37.96	1.50	297.59
321-A [A+B]	140	1.75	0.50	2.00	109.50	3.00	34.00	1.00	296.30
323-A [A+B]	140	1.50	0.50	2.00	109.00	3.00	34.00	1.00	295.55
324-A [A+B]	140	1.25	0.50	2.00	108.50	3.00	34.00	1.00	294.80
326-A [A+B]	139	1.50	0.48	2.50	104.64	2.85	37.96	1.50	297.25
327-A [A+B]	140	1.50	0.50	3.00	110.00	3.00	34.00	1.00	297.55
330-A [A+B]	140	1.75	0.50	3.00	110.50	3.00	34.00	1.00	298.30
331-A [A+B]	140	1.25	0.50	3.00	109.50	3.00	34.00	1.00	296.80
332-A [A+B]	140	1.50	0.50	4.00	111.00	3.00	34.00	1.00	299.55



NIPRO CORPORATION (Headquarters)

3-9-3, Honjo-Nishi, Kita-ku, Osaka,
531-8510, Japan
Tel: +81 6 6372 2331 Fax: +81 6 6371 7422
<https://www.nipro.com.jp/en/>

NIPRO ASIA PTE LTD

6 Temasek Boulevard
#31-01/02 Suntec Tower Four
Singapore 038986
Tel: +65 6588 2281 Fax: +65 6588 2481
Email: niproasia@nipro.com.sg
<http://www.niproasia.com.sg/>

NIPRO MALAYSIA SDN. BHD

B-3-2 The Ascent Paradigm,
No.1 Jalan SS7/26A, Kelana Jaya,
47301 Petaling Jaya, Selangor, Malaysia
Tel: +603 7622 1190/1191/1192
Fax: +603 7622 180
Email: office@nipro.com.my

NIPRO JMI MEDICAL LTD

Unique Heights, Level 11,
117, Kazi Nazrul Islam Avenue, Ramna,
Dhaka-1217, Bangladesh
Tel: +880 2 55138723 Fax: +880 2 55138725
Email: info@jmigroup-bd.com

NIPRO AUSTRALIA PTY LTD

Suite 2. 02, 657 Pacific Highway,
St. Leonards, NSW 2065, Australia
Tel: +61 2 9437 9000 Fax: + 61 2 9437 0860
Email: support@niproaustralia.com.au
<http://www.niproaustralia.com.au/>

NIPRO MEDICAL PVT LTD

24-CCA, DHA Phase 8, (Ex Park View),
Lahore-54792, Pakistan
Tel: +92 4237175802/4237175803 Fax: +92 4237175804
Email: info@nipromedical.com.pk

NIPRO SALES (THAILAND) CO., LTD.

89/170 Nipro Building, Viphavadee-Rangsit Road,
Talad-Bangkhen, Laksi, Bangkok, THAILAND 10210
Tel : +66 2551 2000 Fax : +66 2551 2006
Email : sales@niprosales.com

NIPRO SALES VIETNAM COMPANY LIMITED

Floor 15, Unit 15.06 & 15.07, Viettel Complex Building 285
Cach Mang Thang 8, Ward 12, District 10,
Ho Chi Minh city, Vietnam
Tel : +84 908 260 405
Email : nhuy@niprosales.com

**NIPRO SALES (THAILAND) CO., LTD. -
MYANMAR REPRESENTATIVE OFFICE**

No. E 302, Manawhary Housing, Baho Road, Ahlone
Township,
Yangon, Myanmar
Tel : +66 2551 2000 Fax : +66 2551 2006
Email : win@niprosales.com