



BIOHOLE™ CATHETER NEEDLE

REDUCED PAIN AND INCREASED COMFORT DURING
ARTERIOVENOUS INSERTION



BIOHOLE CATHETER NEEDLE ADVANTAGES

- Plastic fistula needles may prevent vascular access complications 1
- Safe cannulation during fistula maturation 1
- Decreased risk of perforation 1
- Greater patient mobility
- Better option for patients allergic to metal

SIDE HOLES

Optimizes blood flow and minimizes the risk of aspiration to the vessel wall

ULTRA-THIN CATHETER WALL

Guarantees a nearly painless catheter insertion and optimal blood flow

COLOR-CODED HUB

Facilitates the recognition of gauge size

SOFT FLEXIBLE TUBING

Eases clamping during the bloodline connection/disconnection while absorbing shocks and reducing needle movement in patient

TAPERED LUER

Warrants a safe and easy connection/disconnection with bloodlines

SILICONIZED DULL NEEDLE

Ensures smooth penetration and prevents vessel trauma



BIOHOLE NEEDLES

Form the perfect combination with the BioHole Buttonhole Device



SPECIFICATIONS

Catheter Needle*	Size (Gauge)	Catheter Length (mm)	Color
AV-D-CATH-1425-GJ	14	25	Purple
AV-D-CATH-1432-GJ	14	32	Purple
AV-D-CATH-1438-GJ	14	38	Purple
AV-D-CATH-1525-GJ	15	25	Yellow
AV-D-CATH-1532-GJ	15	32	Yellow
AV-D-CATH-1538-GJ	15	38	Yellow
AV-D-CATH-1625-GJ	16	25	Green
AV-D-CATH-1632-GJ	16	32	Green
AV-D-CATH-1638-GJ	16	38	Green
AV-D-CATH-1725-GJ	17	25	Orange
AV-D-CATH-1732-GJ	17	32	Orange
AV-D-CATH-1738-GJ	17	38	Orange

Please contact your local office for other specifications, order quantities, and availability.



NIPRO CORPORATION (Headquarters)

3-9-3, Honjo-Nishi, Kita-ku, Osaka,
531-8510, Japan
Tel: +81 6 6372 2331 Fax: +81 6 6371 7422
<https://www.nipro.com.jp/en/>

NIPRO ASIA PTE LTD

6 Temasek Boulevard
#31-01/02 Suntec Tower Four
Singapore 038986
Tel: +65 6588 2281 Fax: +65 6588 2481
Email: niproasia@nipro.com.sg
<http://www.niproasia.com.sg/>

NIPRO MALAYSIA SDN. BHD

B-3-2 The Ascent Paradigm,
No.1 Jalan SS7/26A, Kelana Jaya,
47301 Petaling Jaya, Selangor, Malaysia
Tel: +603 7622 1190/1191/1192
Fax: +603 7622 180
Email: office@nipro.com.my

NIPRO JMI MEDICAL LTD

Unique Heights, Level 11,
117, Kazi Nazrul Islam Avenue, Ramna,
Dhaka-1217, Bangladesh
Tel: +880 2 55138723 Fax: +880 2 55138725
Email: info@jmigroup-bd.com

NIPRO AUSTRALIA PTY LTD

Suite 2. 02, 657 Pacific Highway,
St. Leonards, NSW 2065, Australia
Tel: +61 2 9437 9000 Fax: + 61 2 9437 0860
Email: support@niproaustralia.com.au
<http://www.niproaustralia.com.au/>

NIPRO MEDICAL PVT LTD

24-CCA, DHA Phase 8, (Ex Park View),
Lahore-54792, Pakistan
Tel: +92 4237175802/4237175803 Fax: +92 4237175804
Email: info@nipromedical.com.pk

NIPRO SALES (THAILAND) CO., LTD.

89/170 Nipro Building, Viphavadee-Rangsit Road,
Talad-Bangkhen, Laksi, Bangkok, THAILAND 10210
Tel : +66 2551 2000 Fax : +66 2551 2006
Email : sales@niprosales.com

NIPRO SALES VIETNAM COMPANY LIMITED

Floor 15, Unit 15.06 & 15.07, Viettel Complex Building 285
Cach Mang Thang 8, Ward 12, District 10,
Ho Chi Minh city, Vietnam
Tel : +84 908 260 405
Email : nhuy@niprosales.com

**NIPRO SALES (THAILAND) CO., LTD. -
MYANMAR REPRESENTATIVE OFFICE**

No. E 302, Manawhary Housing, Baho Road, Ahlone
Township,
Yangon, Myanmar
Tel : +66 2551 2000 Fax : +66 2551 2006
Email : win@niprosales.com