


# NIPROSOL

ACID CONCENTRATE CANISTERS



# NIPROSOL

## Acid concentrate canisters

Liquid acid concentrates intended to be used with powder bicarbonate concentrate.

The ready-to-use dialysis fluid is prepared by the hemodialysis machine by mixing the acid concentrate with the bicarbonate concentrate and the purified water. The standard mixing ratios of the acid concentrates is 1+44.



## NIPROSOL PORTFOLIO

	Na <sup>+</sup>	Ca <sup>++</sup>	Mg <sup>++</sup>	K <sup>+</sup>	Cl <sup>-</sup>	CH <sub>3</sub> COO <sup>-</sup>	HCO <sub>3</sub> <sup>-</sup>	Glucose	Osmolarity
	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	g/l	mOsm/l
306-A [A+B]	139	1.25	0.48	1.50	103.14	2.85	37.96	1.50	294.50
307-A [A+B]	139	1.75	0.48	1.50	104.14	2.85	37.96	1.50	296.00
308-A [A+B]	139	1.50	0.48	1.50	103.64	2.85	37.96	1.00	292.48
310-A [A+B]	139	1.50	0.50	1.00	103.20	2.85	37.96	1.00	291.55
313-A [A+B]	140	1.50	0.50	1.50	108.50	3.00	34.00	1.00	294.55
314-A [A+B]	140	1.25	0.50	1.50	108.00	3.00	34.00	1.00	293.80
316-A [A+B]	140	0.63	0.25	1.50	102.44	4.00	36.81	1.00	291.18
318-A [A+B]	140	1.75	0.50	2.00	109.50	3.00	34.00	1.50	299.08
319-A [A+B]	139	1.75	0.48	3.00	105.64	2.85	37.96	1.50	299.00
320-A [A+B]	139	1.25	0.50	3.00	104.70	2.85	37.96	1.50	297.59
321-A [A+B]	140	1.75	0.50	2.00	109.50	3.00	34.00	1.00	296.30
323-A [A+B]	140	1.50	0.50	2.00	109.00	3.00	34.00	1.00	295.55
324-A [A+B]	140	1.25	0.50	2.00	108.50	3.00	34.00	1.00	294.80
326-A [A+B]	139	1.50	0.48	2.50	104.64	2.85	37.96	1.50	297.25
327-A [A+B]	140	1.50	0.50	3.00	110.00	3.00	34.00	1.00	297.55
330-A [A+B]	140	1.75	0.50	3.00	110.50	3.00	34.00	1.00	298.30
331-A [A+B]	140	1.25	0.50	3.00	109.50	3.00	34.00	1.00	296.80
332-A [A+B]	140	1.50	0.50	4.00	111.00	3.00	34.00	1.00	299.55



**NIPRO CORPORATION (Headquarters)**

3-9-3, Honjo-Nishi, Kita-ku, Osaka,  
531-8510, Japan  
Tel: +81 6 6372 2331 Fax: +81 6 6371 7422  
<https://www.nipro.com.jp/en/>

**NIPRO ASIA PTE LTD**

6 Temasek Boulevard  
#31-01/02 Suntec Tower Four  
Singapore 038986  
Tel: +65 6588 2281 Fax: +65 6588 2481  
Email: [niproasia@nipro.com.sg](mailto:niproasia@nipro.com.sg)  
<http://www.niproasia.com.sg/>

**NIPRO MALAYSIA SDN. BHD**

B-3-2 The Ascent Paradigm,  
No.1 Jalan SS7/26A, Kelana Jaya,  
47301 Petaling Jaya, Selangor, Malaysia  
Tel: +603 7622 1190/1191/1192  
Fax: +603 7622 180  
Email: [office@nipro.com.my](mailto:office@nipro.com.my)

**NIPRO JMI MEDICAL LTD**

Unique Heights, Level 11,  
117, Kazi Nazrul Islam Avenue, Ramna,  
Dhaka-1217, Bangladesh  
Tel: +880 2 55138723 Fax: +880 2 55138725  
Email: [info@jmigroup-bd.com](mailto:info@jmigroup-bd.com)

**NIPRO AUSTRALIA PTY LTD**

Suite 2. 02, 657 Pacific Highway,  
St. Leonards, NSW 2065, Australia  
Tel: +61 2 9437 9000 Fax: + 61 2 9437 0860  
Email: [support@niproaustralia.com.au](mailto:support@niproaustralia.com.au)  
<http://www.niproaustralia.com.au/>

**NIPRO MEDICAL PVT LTD**

24-CCA, DHA Phase 8, (Ex Park View),  
Lahore-54792, Pakistan  
Tel: +92 4237175802/4237175803 Fax: +92 4237175804  
Email: [info@nipromedical.com.pk](mailto:info@nipromedical.com.pk)

**NIPRO SALES (THAILAND) CO., LTD.**

89/170 Nipro Building, Viphavadee-Rangsit Road,  
Talad-Bangkhen, Laksi, Bangkok, THAILAND 10210  
Tel : +66 2551 2000 Fax : +66 2551 2006  
Email : [sales@niprosales.com](mailto:sales@niprosales.com)

**NIPRO SALES VIETNAM COMPANY LIMITED**

Floor 15, Unit 15.06 & 15.07, Viettel Complex Building 285  
Cach Mang Thang 8, Ward 12, District 10,  
Ho Chi Minh city, Vietnam  
Tel : +84 908 260 405  
Email : [nhuy@niprosales.com](mailto:nhuy@niprosales.com)

**NIPRO SALES (THAILAND) CO., LTD. -  
MYANMAR REPRESENTATIVE OFFICE**

No. E 302, Manawhary Housing, Baho Road, Ahlone  
Township,  
Yangon, Myanmar  
Tel : +66 2551 2000 Fax : +66 2551 2006  
Email : [win@niprosales.com](mailto:win@niprosales.com)